




TRAZER for Prosthetists & Physical Therapists

Increasing Efficacy in Fitting & Rehabilitation

**Hosted by:
Randy Cohen, DPT, PT, ATC**



Meet the Panel





Kyle Lorenz
Police Officer |
Tucson Police Department




Jeff Denune, CP/L
Managing Clinical Director |
NuTech Institute

TRAZER

Technology for injury recovery, injury prevention, and enhancing performance, TRAZER tracks, measures, and enhances physical *and* cognitive function.



HIPAA-compliant Data Analytics Portal



100+ Immersive reaction-based activities

4 PILLARS

Activities to Assess and Monitor Essential Metrics

100+ immersive, reaction-based activities in key categories:

Assessments | Rehabilitation | Drills | Workouts | Games



BALANCE

Evaluate stability and adherence to prescribed postures.



DYNAMIC MOVEMENT

Identify and monitor asymmetry in multi-directional movement.



KINEMATICS

Analyze joint range of motion in upper and lower extremities.



NEUROMECHANICS

Improve physical and cognitive accuracy via whole-body decision-making.



BALANCE

BALANCE ASSESSMENT

- How do you determine if balance and stability are affected with the fitting of a prosthetic socket, or a change in ankle, foot, or shoe?
- How do you determine if balance is improving with rehab?
- How do you objectively determine if function is improving?



STEADI Balance


A 4-step Balance Test consisting of progressively challenging tasks designed to assess Users' ability to maintain balance per CDC guidelines.

[Learn More](#)

TRAZER® STEADI BALANCE



STEADI Balance
Modified Tandem



BALANCE

Footwear: Shoe 1
Rigid Foot. Soft Rearfoot and Forefoot.



December 8, 2023

Test 1 | 7:31 AM
Test 2 | 7:32 AM
Test 3 | 7:33 AM


Sway in Feet | Total Movement

- 2.32 | 4 Steps
- 1.59 | 3 Steps
- 3.06 | 1 Step

Subjective: felt stress/pressure points in socket

Sway

Feet (ft)



	Test #1	Test #2	Test #3
Forward	0.77	0.21	1.52
Forward Right	0.29	0.07	0.13
Right	0.00	0.01	0.07
Backward Right	0.00	0.14	0.15
Backward	0.00	0.01	0.49
Backward Left	0.00	0.05	0.16
Left	0.14	0.28	0.31
Forward Left	1.12	0.83	0.22
Total Sway	2.32	1.59	3.06

BALANCE

Sway

Feet (ft)

Footwear: Shoe 2

Solid Rearfoot with side reinforcement with soft Forefoot.




December 8, 2023

Test 1 | 7:44 AM

Test 2 | 7:45 AM

Test 3 | 7:46 AM



Sway in Feet | Total Movement

→ 1.24 | 1 Step

→ 1.18 | 1 Step

→ .79 | 1 Step

Subjective: felt some pressure occasionally in socket

	Test #1	Test #2	Test #3
Forward	0.23	0.00	0.00
Forward Right	0.17	0.00	0.01
Right	0.14	0.00	0.05
Backward Right	0.28	0.54	0.10
Backward	0.23	0.63	0.41
Backward Left	0.14	0.00	0.15
Left	0.02	0.00	0.06
Forward Left	0.03	0.00	0.00
Total Sway	1.24	1.18	0.79

Sway

Feet (ft)

BALANCE

Footwear: Boot

Rigid Foot. Rigid Rearfoot and Forefoot.



December 8, 2023


Test 1 | 7:51 AM

Test 2 | 7:52 AM

Test 3 | 7:53 AM

No misses.

Subjective: felt more stability and less stress on socket



	Test #1	Test #2	Test #3
Forward	0.07	0.22	0.12
Forward Right	0.19	0.09	0.05
Right	0.15	0.02	0.04
Backward Right	0.13	0.00	0.02
Backward	0.00	0.08	0.06
Backward Left	0.00	0.13	0.13
Left	0.00	0.61	0.08
Forward Left	0.00	0.10	0.04
Total Sway	0.55	1.24	0.53


COMPARISON

Best: Shoe 1 Shoe 2 & Shoe 3



Sway

Feet (ft)



	Test #1	Test #2	Test #3
Forward	0.21	0.23	0.12
Forward Right	0.07	0.17	0.05
Right	0.01	0.14	0.04
Backward Right	0.14	0.28	0.02
Backward	0.01	0.23	0.06
Backward Left	0.05	0.14	0.13
Left	0.28	0.02	0.08
Forward Left	0.83	0.03	0.04
Total Sway	1.59	1.24	0.53

- Forward
- Forward Right
- Right
- Backward Right
- Backward
- Backward Left
- Left
- Forward Left
- Total Sway

Misses


	Step Detected	Total
Test #1	3	3
Test #2	1	1
Test #3	0	0

COMPARISON

Best: Shoe 1 vs Shoe 2

Sway

Feet (ft)




	Test #1	Test #2	Difference
Forward	0.21	0.23	+10.68 %
Forward Right	0.07	0.17	+150.49 %
Right	0.01	0.14	+1280.65 %
Backward Right	0.14	0.28	+97.01 %
Backward	0.01	0.23	+2063.64 %
Backward Left	0.05	0.14	+192.47 %
Left	0.28	0.02	-93.81 %
Forward Left	0.83	0.03	-96.37 %
Total Sway	1.59	1.24	-22.07 %

Best: Shoe 1 vs BOOT

Sway

Feet (ft)




	Test #1	Test #2	Difference
Forward	0.21	0.12	-44.90 %
Forward Right	0.07	0.05	-25.98 %
Right	0.01	0.04	+251.61 %
Backward Right	0.14	0.02	-87.59 %
Backward	0.01	0.06	+481.82 %
Backward Left	0.05	0.13	+174.66 %
Left	0.28	0.08	-72.02 %
Forward Left	0.83	0.04	-94.71 %
Total Sway	1.59	0.53	-66.51 %

COMPARISON

Best: **Shoe 2** vs **BOOT**

Sway

Feet (ft)



Forward

Forward Right

Right

Backward Right

Backward

Backward Left

Left

Forward Left

Total Sway

Test #1	Test #2	Difference
0.23	0.12	-50.21 %
0.17	0.05	-70.45 %
0.14	0.04	-74.53 %
0.28	0.02	-93.70 %
0.23	0.06	-73.11 %
0.14	0.13	-6.09 %
0.02	0.08	+351.92 %
0.03	0.04	+45.65 %
1.24	0.53	-57.03 %

0.23

0.12

-50.21 %

0.17

0.05

-70.45 %

0.14

0.04

-74.53 %

0.28

0.02

-93.70 %

0.23

0.06

-73.11 %

0.14

0.13

-6.09 %

0.02

0.08

+351.92 %

0.03

0.04

+45.65 %

1.24

0.53

-57.03 %

CONCLUSION: Shoe Balance Assessment

Measure What Matters

Tested 3 pairs of shoes on same day.

Able to make recommendation of stability based on objective numbers.


Able to get stable forefoot and rearfoot shoes with lower profile than boot with same stability.

DISCUSSION

Impact

- **1 x 30-Minute Session**
- **3 Shoes**
- **80%+ Better Stability**
- **Kyle thoughts**
- **Jeff impression. Use for other components?**

DYNAMIC MOVEMENT

An illustration of a woman in a green t-shirt and black leggings standing on a white square platform. To her right is a circular icon with eight arrows pointing outwards, representing the 8 randomized directions of the drill.

React

A single-step Drill designed to Test or Train reactive agility in 8 randomized directions.

[Learn More](#)



REACT

Movement in 8 Directions from Center

Test 1: June 5, 2023 vs Test 2: December 18, 2023

Step Summary


	Duration	Targets	Calories	Distance	Deceleration Index
Test #1	2:00	12.00	6.68	107.28 ft	1.85
Test #2	6:47	72.00	41.54	754.12 ft	0.85

Averages

	Reaction Time	Dynamic Reaction Time	Speed	Acceleration	Deceleration
Test #1	0.85 s	0.85 s	1.06 ft/s	1.68 ft/s ²	1.68 ft/s ²
Test #2	0.39 s	0.50 s	2.00 ft/s	3.25 ft/s ²	2.93 ft/s ²
Difference	-53.95 %	-42.00 %	+89.99 %	+93.16 %	+74.41 %


Dynamic Reaction Time

Seconds



Speed


ft/s



REACT

Movement in 8 Directions from Center

Distance



Date

Total Distance

06/05/2023

107.28

07/10/2023

135.54

07/10/2023

156.42

08/04/2023

118.91

08/04/2023

120.17

09/15/2023

194.43

11/27/2023

584.95


12/18/2023

754.12

REACT

Movement in 8 Directions from Center

Reaction Time



Date

Avg RT

06/05/2023

0.85

07/10/2023

0.85

07/10/2023

0.48

08/04/2023

0.54

08/04/2023

0.51

09/15/2023

1.31

11/27/2023

1.30


12/18/2023

0.39

REACT

Movement in 8 Directions from Center

Dynamic Reaction Time



Date

Avg Dynamic RT

06/05/2023

0.85

07/10/2023

0.91

07/10/2023

0.51

08/04/2023

0.58

08/04/2023

0.62

09/15/2023

1.36

11/27/2023

1.31


12/18/2023

0.50

REACT

Movement in 8 Directions from Center

Speed



Date

Avg Speed

06/05/2023

1.06

07/10/2023

1.28

07/10/2023

1.43

08/04/2023

1.68

08/04/2023

1.68

09/15/2023

1.81

11/27/2023

1.60


12/18/2023

2.00

REACT

Movement in 8 Directions from Center

Acceleration



Date

Avg Acceleration

06/05/2023

1.68

07/10/2023

1.77

07/10/2023

2.15

08/04/2023

2.94

08/04/2023

2.85

09/15/2023

3.58

11/27/2023

2.56


12/18/2023

3.25

REACT

Movement in 8 Directions from Center

Deceleration



Date

Avg Deceleration

06/05/2023

1.68

07/10/2023

1.74

07/10/2023

2.09

08/04/2023

2.81

08/04/2023

2.72

09/15/2023

3.59

11/27/2023

2.57

12/18/2023

2.93

DISCUSSION


Impact

**Kyle: 1st experience when doing dynamic movement on Trazer?
Experience using when have new devices?**

Jeff: Thoughts on this movement pattern and data?

MINI-T


Measure speed and reactive agility



Mini T

A single-step Drill measuring speed and reactive agility, requires Users to move to random prompts in a Mini-T formation L/R then back to center.

[Learn More](#)



MINI-T

Measure speed and reactive agility

Step Summary

Mini-T	Duration	Targets	Calories
	3:33	12.00	6.48
Distance	Deceleration Index		
111.13 ft	0.36		

Averages

	Reaction Time	Dynamic Reaction Time	Speed	Acceleration	Deceleration
Test #1	0.53 s	0.80 s	2.06 ft/s	3.26 ft/s ²	2.61 ft/s ²


DISCUSSION

Impact

Jeff: Benefits of this movement pattern

Kyle: Unique challenge of this pattern?

KINEMATICS




Squat

0544 -

A 5-rep squat Test to measure lower extremity joint angles.


[Learn More](#)



KINEMATICS

Double-leg Squat

Averages



Test #1

Stance Width Ratio	1.75
Stance Width Distance	1.66 ft
Squat Depth	15.58 in
Trunk Lean	37° F
Pelvic Tilt	5° L
Right Knee Valgus/Varus	8° Varus
Right Knee Flexion	110°
Right Ankle Dorsiflexion	23°
Left Knee Valgus/Varus	11° Varus
Left Knee Flexion	105°
Left Ankle Dorsiflexion	19°

TRAZER O & P ASSESSMENT



O & P Assessment

Assessments

A multi-step Test involving ambulation and balance tasks to measure functional mobility.

1. TUG

Position a chair at the back center of the grid (12.5 feet from the Unit to the center of the seat). Calibrate standing, then sit in the chair over the start point circle. When a bumper appears, stand and walk to the bumper then turn and walk back to the chair and sit. 2 targets

2. Chair Stand-30

Position a chair behind center start point. Calibrate sitting at the edge of the chair with your feet in the start point circle. Stand up when a green ball appears. Sit down when a yellow ring appears. 0:30.

3. Mini-T

From back center start point, move to front center bumper, then move to random left and right bumpers forming a T pattern. Immediately backpedal to the flashing reset point for the next target. 40 targets.

4. Double Leg STEADI Balance

Stand straight with feet shoulder-width apart. Calibrate then hold the stance keeping the on-screen ball as close to the center circle as possible. 0:20.

TRAZER O & P ASSESSMENT



TUG

Based on the CDC Timed Up and Go Test (TUG), this Assessment measures the time it takes the User to stand from a seated position, walk 6 feet, and return to the seated

[Learn More](#)



Chair Stand-30

This single-step Test measures lower body strength and endurance as the User stands up and sits down as many times as possible in 30 seconds.

[Learn More](#)



Mini-T

A Reactive Agility Test measuring speed and agility, requires Users to move to random prompts in a Mini-T formation L/R then back to center.

[Learn More](#)



STEADI Balance

A 4-step Balance Test consisting of progressively challenging tasks designed to assess Users' ability to maintain balance per CDC guidelines.

[Learn More](#)



Q & A




TRAZER for Prosthetists & Physical Therapists

Increasing Efficacy in Fitting & Rehabilitation

**Hosted by:
Randy Cohen, DPT, PT, ATC**



Thank you!



Kyle Lorenz
Police Officer |
Tucson Police Department



Jeff Denune, CP/L
Managing Clinical Director |
NuTech Institute

Thank you!



Randy Cohen
Randy@trazer.com



Jeff Denune
jeff.denune@nutechsyn.com



LEARN MORE
[TRAZER.COM/NUTECH](https://trazer.com/nutech)

SEE US IN CHICAGO AT AAOP

Booth **209**

Technical Workshop
Wednesday, March 6
2:20 PM | Roosevelt 1AB