



SPORTS MEDICINE SUMMIT

Performance Enhancement: From Data to Results

Wednesday, April 17, 2024



MEET THE PANEL



Dr. Randy Cohen

DPT, ATC | PT at Physio Shop | VP at Trazer Inc.

University of Arizona



Dr. Connor Norman

Director of Sports Medicine-Football

University of Georgia





Nathan Wilder

MS, ATC, CSCS | Senior AD, Sports Med - Performance


Towson University

TRAZER

Technology for injury recovery, injury prevention, and enhancing performance, TRAZER tracks, measures, and enhances physical *and* cognitive function.



HIPAA-compliant Data Analytics Portal



100+ Immersive reaction-based activities

All 150+ Activities are built across our 4 Pillars. Each Pillar presents a different experience to the User, analyzes specific types of movement, and reports relevant metrics.

TRAZER: Objective, Comprehensive, and Data-Driven



Balance

WHAT

Balance Activities are built to analyze stability and the ability to maintain and complete a required posture.

HOW

Stand in the required position and stay as still as possible for 20 seconds. The on-screen ball will move when you sway, so try to keep it on the center of the target.

WHY

TRAZER reports how stable you were during the Activity and displays how much you swayed in each direction and the number of times you moved out of position.



Dynamic Movement

WHAT

Dynamic Movement Activities are built to analyze directional movement as the body starts, stops, and changes direction in response to a single-task of moving to the requested target location.

HOW

When you see a target appear in the simulation, move your body to reach the target location, using the avatar to gauge your proximity to the target.

WHY

TRAZER reports how well you move and respond to a changing environment. We will see your reaction time, speed, acceleration, and deceleration and identify any issues in your movement patterns.

Kinematics

WHAT

Kinematics Activities are built to analyze joint angles during repetitions of a specific motion.

HOW

Follow the on-screen instructions to complete 5 repetitions of the movement while TRAZER tracks your joint angles.

WHY

TRAZER reports the angles of your joints so we can see if there are any inconsistencies or imbalances in your range of motion and form.

Neuromechanics

WHAT

Neuromechanics Activities are built to analyze directional movement as the body starts, stops, and changes directions in response to a dual-task of both a cognitive challenge and choice of 2 target locations.

HOW

When you see a prompt appear at the top of the screen, choose between the 2 targets and move to reach its location, using the avatar to gauge your proximity to the target.

WHY

TRAZER reports how well you move in response to the cognitive prompt. We will see your reaction time, speed, acceleration, and deceleration and identify any issues in your movement patterns as well as your cognitive scores.

MAXIMIZING DATA CAPTURE USING TRAZER

Activities to Assess and Monitor Essential Metrics

INTRODUCING THE 4 PILLARS




BALANCE

Direction of Sway |
Errors |
Different levels of
Balance Challenge



DYNAMIC MOVEMENT

Reaction Time | Speed |
Acceleration |
Deceleration | Fatigue
Measure



KINEMATICS

Trunk Lean | Pelvic Drop |
Joint Angle |
Valgus/Varus



NEUROMECHANICS

Reaction Time | Dynamic
Reaction Time | Cognitive
Assessment |
Acceleration |
Deceleration | Fatigue
Measure

BASELINING USING TRAZER

4 PILLARS COMPREHENSIVE BASELINES

Dynamic Movement

Balance

Kinematics

Neuromechanics



Bel L, Mathieu N, Ducrest V, Bizzini M. Lower Limb Exercise-Based Injury Prevention Programs Are Effective in Improving Sprint Speed, Jumping, Agility and Balance: an Umbrella Review. *IJSPT*. 2021;16(6):1396-1404. [doi:10.26603/001c.29860](https://doi.org/10.26603/001c.29860)

Umbrella Review

Lower Limb Exercise-Based Injury Prevention Programs Are Effective in Improving Sprint Speed, Jumping, Agility and Balance: an Umbrella Review

Loïc Bel, PT, Bsc Msc Cand.^{1, a}, Nicolas Mathieu, PT, BSc, MSc¹, Vincent Ducrest, PT, BSc, MSc¹

¹ Physiotherapy, HES-SO Valais, University of Applied Sciences Western Switzerland, ² Research, Schulthes Institute, University of Applied Sciences Western Switzerland

Keywords: systematic review, performance, prevention programs, lower limb injury, exercise

<https://doi.org/10.26603/001c.29860>



Journal of Sport and Health Science

Volume 1, Issue 1, May 2012, Pages 36-42



International Journal of Sports Physical Therapy

Vol. 16, Issue 6, 2021

16(6):1396-1404

International Journal of Sports Physical Therapy

<https://doi.org/10.26603/001c.29860>

16(6):1396-1404

ORIGINAL RESEARCH article

Front. Sports Act. Living. 01 October 2021
Sec. Injury Prevention and Rehabilitation
Volume 3 - 2021 | <https://doi.org/10.3389/fspor.2021.729729>

A Novel Approach to Assessment of Perceptual-Motor Efficiency and Training-Induced Improvement in the Performance Capabilities of Elite Athletes



Gary B. Wilkerson^{1*}



Dustin C. Nabhan^{2*}



Tyler S. Perry^{3*}

Original article

Biomechanical risk factors of non-contact ACL injuries: A stochastic biomechanical modeling study

Cheng-Feng Lin^a, Hui Liu^b, Michael T. Gros^c, Paul Weinholt^{d,e}, William E. Garrett^f, Bing Yu^{c,d,e}  

Assessment and Training of Perceptual-Motor Efficiency for Human Performance Optimization

[View all 3 Articles >](#)

Impact of Baseline Testing on Concussion Recovery

4 STEPS **TO BASELINES**

Cognitive Function
Assessment

Balance Assessment

Physical Health Evaluation

Objective Measurements

Baseline testing: 11.7 days recovery
(vs. 19.8 days without)

Improved recovery time = 8.1 days

Clinical Journal of Sport Medicine

Baseline testing: 29.3% prolonged recovery
(vs. 42.9% without)

Improved recovery rate = 13.6%

Journal of Athletic Training

**Nathan Wilder - Senior Athletic Director,
Towson University**



Nathan Wilder

**MS, ATC, CSCS | Senior AD,
Sports Med - Performance**

Towson University


TOWSON BASELINE



BESS Test

A multi-step balance Test of Bilateral, Unilateral and Tandem postures conducted on firm and foam surfaces.


[Learn More](#)



Linear

A single-step Test where the User moves forward to a target for 10 repetitions to measure the speed and efficiency of movement.

[Learn More](#)



Lateral

A single-step Test that measures speed, acceleration/deceleration and L/R asymmetries as the User moves as fast as possible laterally between pre-defined targets.

[Learn More](#)



Squat

A 5-rep squat Test to measure lower extremity joint angles.

[Learn More](#)



LAS 20

20 randomly appearing targets to the Left and Right of the User allow for isolated measurement of dynamic reaction time and movement symmetry.

[Learn More](#)



Stroop

A single-step dual-task Test designed to measure reaction time in relation to cognitive function. When a word in a certain color appears on screen, the User is required to select between

[Learn More](#)



Single Leg Squat

A multi-step single-leg squat Test to compare lower extremity joint angles of the left and right legs for 5 reps each.

[Learn More](#)


Towson: Pre-Participation Exam



Balance 1



Balance 2



Kinematics



Injury to the First MTP Joint

LAS-20







Linear



Lateral

- **Presentation:** Wide Receiver pushed off L foot while making a block and felt pain in L 1st MTP joint
- **MRI Findings:**
 - Grade 3 sprain of medial collateral ligament
 - Grade 2 sprain of distal plantar plate at 1st MTP joint
 - Grade 2 flexor hallucis muscle strain
- **X-ray Results:** X-rays did not show any sesamoid retraction so opted for conservative management
- **Treatment Plan:**
 - 8 weeks in boot
 - 4-week progression back to Activities of Daily Living (ADLs) in shoes with orthotics
 - 4-week functional progression
 - 4-week progression back to full Return to Play (RTP)





Lateral Agility Screen (LAS) Data

Ability to objectively track week to week improvement of

- Acceleration
- Deceleration
- Speed
- Distance
- Track improvements in overall speed/endurance throughout the test


Deceleration



Speed




Speed



Sagittal Agility Screen Data


Acceleration



Distance




Linebacker Concussion Management




- **Incident:** Linebacker sustained a concussion in a game
- **Assessment and Evaluation:**
 - Removed from play and put through SCAT 5 evaluation by team physician
 - Initial reported 18/22 symptoms with symptom severity score of 38/132
- **Initial Treatment:**
 - Started daily Balance Error Scoring System (BESS) and Stroop tests on TRAZER
 - Added low-level cardio on bike below symptom threshold
- **Monitoring and Tracking:**
 - Tracked TRAZER data and daily symptoms until symptom-free
- **Return to Play (RTP) Protocol:**
 - Initiated concussion RTP progression once symptom-free




BESS Data



Bilateral




Single Leg



Tandem Stance


Dynamic Reaction Time




Dynamic Reaction Time



Speed



Acceleration



Deceleration



Distance



STROOP Data

Grade 3 High Ankle Sprain




- **Incident:** Wide Receiver tackled with a hip drop tackle, resulting in linebacker landing directly onto athlete's foot, forcing dorsiflexion and eversion mechanism of injury (MOI)
- **MRI Findings:**
 - Complete tear of anterior tibiofibular ligament and posterior tibiofibular ligament
 - Distal sprain of the distal interosseous membrane
- **Treatment:**
 - Surgical repair with tightrope procedure and syndesmotic screw fixation
- **Rehabilitation Protocol:**
 - 6 weeks non-weight bearing
 - 4 weeks weight bearing in boot
 - Begin motion and progressive activities as tolerated




Single Leg Squat Kinematic Assessment


Dorsiflexion




Knee Valgus/Varus



Torso




Knee Flexion



Single Leg Squat Kinematic Assessment

Averages




Test #1

Trunk Lean	39° F
Pelvic Tilt	1° R
Right Knee Valgus/Varus	4° Varus
Right Knee Flexion	63°
Right Ankle Dorsiflexion	34°
Left Knee Valgus/Varus	5° Varus
Left Knee Flexion	55°
Left Ankle Dorsiflexion	27°

Dr. Connor Norman
UGA Sports Medicine Director

GEORGIA[®]
ATHLETICS



Dr. Connor Norman

**Director of Sports
Medicine-Football**

University of Georgia

University of Georgia

Objective Measures for Intermediary Assessment

- Utilized **LAS 20** for intermediary assessment before transitioning to field work
- Included objective measures such as:
 - **Squat**
 - **Single Leg (SL) Activities**
 - **Deadlift (DL)**
 - **Jumping, with particular emphasis on Single Leg (SL) jump**
 - Objective measure focused on **trunk lean**
- Endurance metrics tracked with objective numbers
- **Linear and Lateral 2-minute drills** utilized to assess progress over time
- Integrated Catapult data with TRAZER
 - Catapult provided overall acceleration/deceleration metrics
 - TRAZER identified asymmetries in different movement directions



UGA: LAS 20 Report

TRAZER

Step Summary


LAS 20	Duration 1:09	Targets 20.00	Calories 7.04
Distance 129.25 ft	Deceleration Index 1.53		

Averages

	Reaction Time	Dynamic Reaction Time	Speed	Acceleration	Deceleration
Test #1	0.97 s	1.14 s	2.74 ft/s	8.39 ft/s ²	7.92 ft/s ²


TRAZER

Speed



Test #1
Right: 2.84
Left: 2.65
L/R Difference: 6.62% L

Acceleration



Test #1
Right: 8.72
Left: 8.14
L/R Difference: 6.60% L

Deceleration



Test #1
Right: 7.88
Left: 7.95
L/R Difference: 0.92% R

Reaction Time



Test #1
Right: 1.00
Left: 0.93
L/R Difference: 6.66% R

Dynamic Reaction Time



Test #1
Right: 1.18
Left: 1.10
L/R Difference: 6.85% R

Test #1

	Rep #1	Rep #2	Rep #3	Rep #4	Rep #5	Rep #6	Rep #7	Rep #8	Rep #9	Rep #10	Rep #11	Rep #12	Rep #13	Rep #14	Rep #15	Rep #16	Rep #17	Rep #18	Rep #19	Rep #20
Avg Reaction Time	0.94 s	0.97 s	1.08 s	0.97 s	0.91 s	0.95 s	0.95 s	1.04 s	1.04 s	1.01 s	0.98 s	0.91 s	1.04 s	1.01 s	0.87 s	1.02 s	0.90 s	0.98 s	0.87 s	0.91 s
Avg Dynamic Reaction Time	0.94 s	0.97 s	1.08 s	0.97 s	0.95 s	2.71 s	0.95 s	1.04 s	1.04 s	1.01 s	0.98 s	0.91 s	1.04 s	1.01 s	0.87 s	1.02 s	0.90 s	0.98 s	2.45 s	0.91 s
Avg Speed	3.77 ft/s	2.47 ft/s	2.75 ft/s	2.33 ft/s	3.32 ft/s	3.36 ft/s	2.48 ft/s	2.33 ft/s	2.79 ft/s	2.95 ft/s	2.40 ft/s	2.44 ft/s	2.69 ft/s	2.60 ft/s	2.56 ft/s	2.88 ft/s	2.91 ft/s	2.75 ft/s	2.87 ft/s	3.42 ft/s
Avg Acceleration	9.88 ft/s²	6.98 ft/s²	7.92 ft/s²	6.44 ft/s²	9.68 ft/s²	10.01 ft/s²	6.62 ft/s²	8.23 ft/s²	13.74 ft/s²	10.23 ft/s²	6.75 ft/s²	6.55 ft/s²	9.73 ft/s²	7.16 ft/s²	6.95 ft/s²	11.03 ft/s²	8.73 ft/s²	8.49 ft/s²	8.43 ft/s²	9.14 ft/s²
Avg Deceleration	9.26 ft/s²	7.86 ft/s²	10.69 ft/s²	7.57 ft/s²	9.84 ft/s²	10.54 ft/s²	9.39 ft/s²	6.16 ft/s²	6.42 ft/s²	18.77 ft/s²	6.50 ft/s²	7.33 ft/s²	6.80 ft/s²	7.06 ft/s²	8.46 ft/s²	6.70 ft/s²	8.98 ft/s²	6.57 ft/s²	9.04 ft/s²	7.74 ft/s²
Total Distance	6.55 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	11.53 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	6.00 ft	10.04 ft	6.00 ft



UGA Athletes on the MOVE





Dr. Randy Cohen

TRAZER VP of Clinical Education and Elite Performance



**Randy P. Cohen AT, PT, DPT
Retired Assoc AD Medical Services
University of Arizona**



The PhysioShop AT/PT/Wellness Center

DYNAMIC MOVEMENT-REACT

Single-step **Drill** in TRAZER

Objective: Enhance Reactive
Agility


How: Test and Train in 8
Directions



REACT

Movement in 8 Directions from Center

Distance



Date

Total Distance

06/05/2023

107.28

07/10/2023

135.54

07/10/2023

156.42

08/04/2023

118.91

08/04/2023

120.17

09/15/2023

194.43

11/27/2023

584.95


12/18/2023

754.12

REACT

Movement in 8 Directions from Center

Reaction Time



Date

Avg RT

06/05/2023

0.85

07/10/2023

0.85

07/10/2023

0.48

08/04/2023

0.54

08/04/2023

0.51

09/15/2023

1.31

11/27/2023

1.30


12/18/2023

0.39

REACT

Movement in 8 Directions from Center

Dynamic Reaction Time



Date

Avg Dynamic RT

06/05/2023

0.85

07/10/2023

0.91

07/10/2023

0.51

08/04/2023

0.58

08/04/2023

0.62

09/15/2023

1.36

11/27/2023

1.31


12/18/2023

0.50

REACT

Movement in 8 Directions from Center

Speed



Date

Avg Speed

06/05/2023

1.06

07/10/2023

1.28

07/10/2023

1.43

08/04/2023

1.68

08/04/2023

1.68

09/15/2023

1.81

11/27/2023

1.60

12/18/2023

2.00

REACT

Movement in 8 Directions from Center

Acceleration



Date

Avg Acceleration

06/05/2023

1.68

07/10/2023

1.77

07/10/2023

2.15

08/04/2023

2.94

08/04/2023

2.85

09/15/2023

3.58

11/27/2023

2.56


12/18/2023

3.25

REACT

Movement in 8 Directions from Center

Deceleration



Date

Avg Deceleration

06/05/2023

1.68

07/10/2023

1.74

07/10/2023

2.09

08/04/2023

2.81

08/04/2023

2.72

09/15/2023

3.59

11/27/2023



2.57

12/18/2023

2.93

KINEMATICS- DOUBLE LEG SQUAT


- Measure Lower Extremity Joint Angles
- 5-Rep Squat Test
- Comprehensive Lower Body Assessment



KINEMATICS

Double-leg Squat

Averages




Test #1

Stance Width Ratio	1.75
Stance Width Distance	1.66 ft
Squat Depth	15.58 in
Trunk Lean	37° F
Pelvic Tilt	5° L
Right Knee Valgus/Varus	8° Varus
Right Knee Flexion	110°
Right Ankle Dorsiflexion	23°
Left Knee Valgus/Varus	11° Varus
Left Knee Flexion	105°
Left Ankle Dorsiflexion	19°

KINEMATICS

MEASURE DOUBLE LEG SQUAT LOADED



Stance Width Ratio

Stance Width Distance

Trunk Lean

Pelvic Tilt

Right Knee Valgus/Varus

Right Knee Flexion

Right Ankle Dorsiflexion

Left Knee Valgus/Varus

Left Knee Flexion

Left Ankle Dorsiflexion

	Test #1	Test #2
Stance Width Ratio	1.30	1.07
Stance Width Distance	1.11 ft	1.07 ft
Trunk Lean	42° F	27° F
Pelvic Tilt	0° R	0° R
Right Knee Valgus/Varus	3° Valgus	8° Valgus
Right Knee Flexion	93°	72°
Right Ankle Dorsiflexion	31°	28°
Left Knee Valgus/Varus	11° Varus	3° Varus
Left Knee Flexion	93°	72°
Left Ankle Dorsiflexion	34°	34°



FUTURE RELEASES FOR XP CUSTOMERS

Double Leg Jump

A **5-rep jump Test** to measure lower extremity joint angles during landing.


Single Leg Jump

A **multi-step single leg jump Test** to compare lower extremity joint angles during landing of the left and right legs for 5 reps each.



- Enhanced Biomechanical Evaluation
- **Asymmetry Detection**
- Comprehensive Performance Monitoring

Double Leg Jump




Pre-Jump
 Post-Jump

Stance Width Ratio
Stance Width Distance
Squat Depth
Jump Height
Trunk Lean
Pelvic Tilt
Right Knee Valgus/Varus
Right Knee Flexion
Right Ankle Dorsiflexion
Left Knee Valgus/Varus
Left Knee Flexion
Left Ankle Dorsiflexion

Test #1

0.97
0.89 ft
12.58 in
12.55 in
42° F
3° L
5° Varus
98°
28°
7° Varus
98°
28°



Pre-Jump
 Post-Jump

Stance Width Ratio
Stance Width Distance
Squat Depth
Jump Height
Trunk Lean
Pelvic Tilt
Right Knee Valgus/Varus
Right Knee Flexion
Right Ankle Dorsiflexion
Left Knee Valgus/Varus
Left Knee Flexion
Left Ankle Dorsiflexion

Test #1


0.97
0.92 ft
11.62 in
12.55 in
27° F
2° L
6° Varus
91°
24°
8° Varus
93°
23°




Reps

Test #1

	Rep #1	Rep #2	Rep #3	Rep #4	Rep #5
Stance Width Ratio	1.10	1.29	1.34	1.23	0.97
Pre Stance Width Distance	1.38 ft	1.43 ft	1.61 ft	1.93 ft	1.66 ft
Post Stance Width Distance	1.21 ft	1.41 ft	1.52 ft	1.39 ft	1.15 ft
Pre Squat Depth	14.67 in	15.08 in	15.34 in	14.17 in	11.71 in
Post Squat Depth	2.16 in	2.11 in	3.23 in	2.92 in	1.49 in
Jump Height	19.02 in	18.21 in	19.15 in	15.39 in	10.83 in
Pre Trunk Lean	43° F	46° F	48° F	53° F	46° F
Post Trunk Lean	9° F	5° F	18° F	18° F	13° F
Pre Pelvic Tilt	10° L	5° L	7° L	9° L	8° L
Post Pelvic Tilt	9° R	5° L	5° L	3° L	6° L
Pre Right Knee Valgus/Varus	4° Varus	9° Varus	8° Varus	8° Varus	8° Varus
Post Right Knee Valgus/Varus	0° Varus	0° Varus	2° Varus	1° Varus	0° Varus
Pre Right Knee Flexion	75°	111°	114°	77°	94°
Post Right Knee Flexion	50°	52°	60°	56°	46°
Pre Right Ankle Dorsiflexion	33°	30°	33°	46°	24°
Post Right Ankle Dorsiflexion	20°	18°	23°	22°	3°
Pre Left Knee Valgus/Varus	5° Varus	8° Varus	12° Varus	10° Varus	11° Varus
Post Left Knee Valgus/Varus	1° Varus	3° Varus	4° Varus	4° Varus	2° Varus
Pre Left Knee Flexion	87°	111°	111°	62°	98°
Post Left Knee Flexion	47°	55°	55°	58°	42°
Pre Left Ankle Dorsiflexion	26°	31°	31°	42°	26°
Post Left Ankle Dorsiflexion	24°	25°	24°	23°	6°



FATIGUE MEASURE?



CONSISTENT LEAN AWAY


HIGHER LEVEL TRAINING & ANALYSIS OF DATA



Lateral

A single-step Test that measures speed, acceleration/deceleration and L/R asymmetries as the User moves as fast as possible laterally between pre-defined targets.

[Learn More](#)



Flanker 2

Assessments

A single-step dual-task Test designed to measure reaction time in relation to cognitive function. 5 arrows appear in multiple locations on the screen pointing randomly to the left or right. The User must move to the target in the direction indicated by the center arrow.



Drop-Step

A single-step Drill Tests or Trains reactive-based change of direction capability to back diagonal targets.

[Learn More](#)



Math Add Subtract Reach

A single-step Drill designed to Test or Train dual-task reaction time where the User must solve math addition and subtraction problems by choosing between two options and

[Learn More](#)



Jump Explosion

The User engages in lateral movements and jumps to return falling balls to a moving conveyor belt. This Game challenges the User's timing and strategy with 7 levels of

[Learn More](#)

Lateral-Endurance as a Fatigue Measure

Measures:

- Speed
- Acceleration/Deceleration
- Left/Right Asymmetries

User moves laterally between pre-defined targets.

Value:

Assess and enhance lateral agility, speed, and symmetry.
Gain insights into movement patterns for improvement.




Date 07/31/2023, 1:14 PM
 Duration 2:14
 Location The Physio Shop
 Status Complete
 Scripted Activity Lateral





Provider Randy Cohen
 Date 09/20/2023, 5:14 PM
 Duration 2:14
 Location The Physio Shop
 Status Complete
 Scripted Activity Lateral

INTERPRETING THE REPORT

Timeline




Timeline



NEUROMECHANICS

Flanker 2



Flanker 2

Assessments

A single-step dual-task Test designed to measure reaction time in relation to cognitive function. 5 arrows appear in multiple locations on the screen pointing randomly to the left or right. The User must move to the target in the direction indicated by the center arrow.



FLANKER 2

WITH EXTERNAL DEVICES



Flanker 2-Report Summary

Step Summary

Flanker 2 95.00% Cognitive	Duration 2:33	Targets 40.00	Calories 14.99	Distance 339.39 ft	Deceleration Index 1.23
--------------------------------------	-------------------------	-------------------------	--------------------------	------------------------------	-----------------------------------

Cognitive

Test #1

Prompt Time
200 ms

Response Accuracy

	Correct	Incorrect
Right Congruent	10	0
Right Incongruent	10	0
Left Congruent	9	1
Left Incongruent	9	1
Total	38	2

Response Time (s)

	Correct	Incorrect
Right Congruent	2.38	2.56
Right Incongruent	2.26	2.76
Left Congruent	2.10	--
Left Incongruent	2.58	--

Reaction Time (s)

	Correct	Incorrect
Right Congruent	0.35	0.16
Right Incongruent	0.19	0.12
Left Congruent	0.59	--
Left Incongruent	0.20	--

Dynamic Reaction Time (s)



	Correct	Incorrect
Right Congruent	1.28	0.16
Right Incongruent	1.06	1.96
Left Congruent	1.25	--
Left Incongruent	1.31	--

Averages


	Reaction Time	Dynamic Reaction Time	Speed	Acceleration	Deceleration
Test #1	0.32 s	1.14 s	2.76 ft/s	5.92 ft/s ²	5.70 ft/s ²

Advanced Balance

Math Reach



TRAZER GAME: Jump Explosion



Objective: Jump Explosion is a game designed to engage users in lateral movements and jumps.

Value: This game challenges the user's timing and strategy with 7 levels of increasing difficulty. It enhances agility, coordination, and reaction time.

How: Balls descend from a conveyor belt in a chain. The user moves under the chain of balls and jumps to push them back up to the conveyor belt before they hit the ground.



Q & A

CONTACT

Questions About Clinical & Performance Applications

Randy Cohen:

Randy@trazer.com

Sales

Craig Waters

Cwaters@trazer.com

www.trazer.com



Dr. Randy Cohen

DPT, ATC | PT at Physio Shop | VP at Trazer Inc.

University of Arizona